


FILETS TUBULAIRES Réf. FMP / FMS / FMIS / FTC / FTH

Description:	Elastic net bandage
--------------	---------------------

Final dimensions:	25 m x2,4 cm	
Composition:	Polyamide,Polyurethane.	
Size variation:	± 10 %	
Threads sq/cm:	421-004	
Latex:	Free	
PVC:	Free	
Re use:	No	
Sterilizable by:	EO	
Colour:	white	
CFU/g	<10^2	
Fluorescence	conform	
Size	4	



Weight (gr) of packing material used in 100 pcs

CE Class: I CE mark: CE Directive: MDD 93/42/EEC amended by 2007/47/EC Applicable norms: EN ISO 13485:2012; EN ISO 9001:2008; ISO 15223-1:2012; EN 1041+A1:2013; EN ISO 11737-1:2006; EN ISO 14971:2012; EN ISO 10993-1/AC:2010; BP XX H-F:2012 EN 980:2008; EN 1041:2008; ENISO11737-1:2006/AC:2009;	Nr and Revision Date:	6 2013/12/05	
Directive: MDD 93/42/EEC amended by 2007/47/EC Applicable norms: EN ISO 13485:2012 ; EN ISO 9001:2008 ; ISO 15223-1:2012 ; EN 1041+A1:2013 ; EN ISO 11737-1:2006 ; EN ISO 14971:2012 ; EN ISO 10993-1/AC:2010 ; BP XX H-F:2012	CE Class:		
Applicable norms: EN ISO 13485:2012 ; EN ISO 9001:2008 ; ISO 15223-1:2012 ; EN 1041+A1:2013 ; EN ISO 11737-1:2006 ; EN ISO 14971:2012 ; EN ISO 10993-1/AC:2010 ; BP XX H-F:2012	CE mark:	CE	
EN ISO 11737-1:2006; EN ISO 14971:2012; EN ISO 10993-1/AC:2010; BP XX H-F:2012	Directive:	MDD 93/42/EEC amended by 2007/47/EC	
	Applicable norms:	EN ISO 13485:2012; EN ISO 9001:2008; ISO 15223-1:2012; EN 1041+A1:2013;	
EN 980:2008 : EN 1041:2008 : ENISO11737-1:2006/AC:2009 ;		EN ISO 11737-1:2006; EN ISO 14971:2012; EN ISO 10993-1/AC:2010; BP XX H-F:2012;	

PACKING		With IFU- Instructions for use			
Level	Units	Туре	Dimension	Stamp	Labels
1	1	Dispenser shelf box		stamped	labeled
2	8	Carton box	270x145x248mm - BV	stamped	labeled

Important Notes: Nr pcs/pack subjected to ± 3% variation	
	Store in a dry place at room temperature and protect from sunlight. The product has a good stability and
Storage instructions therefore normal conditions of medical devices storage shall apply	

1			
	Waste instructions	medical waste.	
		No special requirements. The waste must be disposed according to local legal regulations for the disposable of	
		therefore,normal conditions of medical devices storage shall apply	
Ш		Store in a dry place at room temperature and protect from sunlight. The product has a good stability and	

ISSUED BY: Quality Management Department, 2015/08/28

WASTE AND ENVIRONMENT REGULATIONS

CONFIDENTIAL: The information contained herein is confidential and belongs to Bastos Viegas, s.a. It is provided for specific purposes and in connection with established agreements. It is not permitted disclosure, reproduction, copying or distribution, in whole or in part, for other purposes or to any third parties.



11, avenue Marcel Dassault – F 37200 TOURS TECHNOPOLE Tél. : 33 (0) 247 28 29 30 – Fax. : 33 (0) 247 28 30 40





5308